

**COURSE TITLE : PHYSICAL CHEMISTRY
PRACTICALS**

COURSE CODE : 15U6PRCHE06

**UNIT 1 : PHYSICAL CHEMISTRY
PRACTICALS**

SESSION 2 : Conductometric Titrations

CONDUCTOMETRIC TITRATIONS

- Conductometric titration is quantitative analysis used to estimate the concentration of a given analyte solution.
- Conductometric titration involves the continuous addition of a titrant to an analyte solution and measuring the conductance of the solution.
- The conductance of the solution depends on the number of free ions and the mobility of ions.
- The conductance of a solution can be measured with the help of a conductivity meter and conductivity cell.
- The end point in a conductometric titration is obtained by plotting a **conductance vs volume of titrant** graph.

STRONG ACID Vs STRONG BASE TITRATION

AIM:

To determine the strength of a given solution of HCl by conductometric titration with the given standard NaOH solution.

PRINCIPLE

Reaction between NaOH and HCl:



- Conductance of acid soln* – Prominently due to H^+ ions
- Addition of alkali – Removes H^+ ions as undissociated H_2O
- Conductance decreases
- After equivalence point - Conductance increases due to free OH^- ions

MATERIALS REQUIRED:

1. Conductivity Meter and Conductivity Cell





2. 0.1 M NaOH
3. Unknown HCl
4. Beaker – 250 mL
5. Burette
6. Pipette – 10 mL
7. Standard Flasks – 100 mL (2 Nos)

PROCEDURE:

- The titrant (Standard NaOH solution) is filled in the burette.
- Pipette out 10 mL of unknown HCl solution into a 250 mL beaker and dilute to 50 mL.
- Dip the conductivity cell in the solution and note the conductance of the solution.
- Add 1 mL of standard NaOH solution from the burette, stir the solution and note the conductance.
- The addition of titrant is repeated in amounts of 1 mL and the conductance is measured each time.

GRAPH

Conductance (*y-axis*) and Volume of titrant added (*x-axis*).



CALCULATION

Volume of unknown HCl (V_1) = 10 mL

Molarity of NaOH (M_2) = 0.1 M

Volume of NaOH (V_2) = ----- mL (determined from graph)

$$\text{Molarity of HCl } (M_1) = \frac{M_2 \times V_2}{V_1} = \frac{0.1 \times V_2}{10}$$