

**COURSE TITLE : PHYSICAL CHEMISTRY
PRACTICALS**

COURSE CODE : 15U6PRCHE06

SESSION 2 : Potentiometric Titrations

POTENTIOMETRIC TITRATIONS

- Potentiometric titration is quantitative analysis used to estimate the concentration of a given analyte solution.
- Potentiometric titration involves the continuous addition of a titrant to an analyte solution and measuring the EMF of the cell.
- The equivalence point of the titration reaction will be indicated by the sudden change in the potential.
- The potential of an electrode varies with ion concentration in accordance with

Nernst equation:

$$E = E^0 + \frac{2.303RT}{nF} \log C$$

E is the electrode potential at the concentration C


KMnO₄ Vs Fe²⁺ soln

AIM:

To determine the strength of a given ferrous ion solution by potentiometric titration with the given KMnO₄ solution.

PRINCIPLE

The redox reaction for the KMnO_4 vs Fe^{2+} titration:



A platinum electrode - ferrous salt solution - calomel electrode

Calomel || KCl || Fe^{3+} , Fe^{2+} | Pt



Indicator electrode – Platinum electrode



Reference electrode – Calomel electrode


The potential acquired by the indicator electrode:

$$E = E^0 + \frac{0.0591}{n} \ln \frac{[Fe^{3+}]}{[Fe^{2+}]}$$

- The EMF of the cell depends on the ratio $[Fe^{3+}]/[Fe^{2+}]$.
- In the beginning, the potential depends on the concentration of Fe^{2+} state.
- After the equivalence point, the potential depends on the concentration of Fe^{3+} state.
- The potential shows a sharp increase at the equivalence point.

MATERIALS REQUIRED:

1. Potentiometer, Platinum and Calomel electrodes



2. 0.1 N KMnO_4

3. Unknown Ferrous ion solution

4. 4N Sulphuric acid

5. Burette

6. Pipette – 10 mL

7. Standard Flasks – 100 mL (2 Nos)

PROCEDURE:


- The unknown Ferrous ion solution is made upto 100 mL
- The titrant (0.1 N KMnO_4 solution) is filled in the burette.
- Pipette out 20 mL of unknown Ferrous ion solution into a 250 mL beaker.
- Add 40 mL of 6N sulphuric acid to the ferrous ions solution.
- Dip the platinum and calomel electrodes in the solution and connected to the potentiometer.
- Record the EMF of the solution.
- Add 1 mL of 0.1 N KMnO_4 solution from the burette, stir the solution and note the EMF.
- The addition of titrant is repeated in amounts of 1 mL and the EMF is measured each time.

$$\frac{\Delta E}{\Delta V}$$

Volume of KMnO_4 (V) added (in mL)	EMF (E) (in Volts)	ΔE	$\frac{\Delta E}{\Delta V}$
0	330	-	-
1	350	20	20
2	360	10	10
3	380	20	20
4			
5			
6			
7			
8			
9			

GRAPH

$\frac{\Delta E}{\Delta V}$ (y-axis) and Volume of titrant added (x-axis).



CALCULATION

Volume of unknown Ferrous ion soln (V_1) = 20 mL

Strength of KMnO_4 soln (N_2) = 0.1 N

Volume of unknown Ferrous ion soln (V_2) = ----- mL (determined from graph)

Normality of Fe^{2+} solution (N_1) = $\frac{N_2 \times V_2}{V_1} = \frac{0.1 \times V_2}{20}$

RESULT:

The strength of the given Fe^{2+} solution isN.