

The background features a smooth gradient from bright yellow at the top to a vibrant green at the bottom. Scattered across this background are numerous water droplets of varying sizes. Some droplets are large and prominent, while others are small and subtle. Each droplet has a realistic 3D effect with a bright highlight on its upper left surface and a soft shadow on its lower right, giving them a sense of depth and volume.

TISSUE CULTURE PROTOCOL

The background features a smooth gradient from bright yellow at the top to a vibrant green at the bottom. Scattered across this background are numerous water droplets of varying sizes. Some droplets are large and prominent, while others are small and subtle. Each droplet has a realistic highlight and shadow, giving them a three-dimensional appearance.

A. TISSUE CULTURE LABORATORY

1. A WASHING AREA

- Large sink
- Running hot and cold tap water
- Brushes of various sizes
- Detergent
- Single distilled water
- Plastic buckets
- Plastic bucket with lid

2. HOT AIR OVEN

- Necessary for drying the washed glass goods.
- Number of enameled trays of different sizes to keep glass goods inside the oven



3. REFRIGERATOR

- For storing various chemical like vitamins, hormones, amino acids, yeast extract, coconut milk, casein-hydrolysate, Stock solutions of salts, etc.



4. DISTILLATION PLANT

- Single distillation and a double distillation water plant.
- Plastic containers for storing the distilled water



5. WEIGHING BALANCE

- Pan balance
- Chemical balance
- Electronic balance

- For weighing chemicals, sugars, agar-agar, etc.



6. pH METER

- For adjusting pH of the nutrient medium



7. AUTOCLAVE

- Important for sterilization of nutrient media, glass goods, instruments, etc.



8. Other Requirements

- **Vacuum Pump**
- **Working Tables**
- **Heater**
- **Microscopes – simple, compound, inverted binocular, camera attachment for photomicrograph**
- **Microtome**
- **Wooden Racks**

9. LABORATORY FOR ASEPTIC INOCULATION

- No windows, no ventilator to make it dust-free.
- Double doors
- Automatic door closer
- Inner door - Rubber mat for cleaning
- Simple inoculating hood

9. LABORATORY FOR ASEPTIC INOCULATION

LAMINAR AIR FLOW

- Number of small blower motors
- Air passes through number of HEPA (High Efficiency Particulate Air) filters
- $0.3\mu\text{m}$
- Ultra clean air free from fungal/ bacterial contaminants
- Flows at a velocity of $27\pm 3\text{m/minute}$
- Put on for 10-15 min before starting work



10. CULTURE ROOM

- For incubating the culture under controlled temperature, light and humidity
- Double doors
- Dust free
- Constant room temperature $25 \pm 2^{\circ}\text{C}$
- Shelves made of glass, plywood
- LED lamps with timer
- Thermometer and Hygrometer (relative humidity should be above 50%)
- Shaker for suspension culture





B. GLASS GOODS AND INSTRUMENTS

1. Glass goods

- Should be of Corning or Pyrex or similar Boro-silicate glass
- Measuring cylinders
- Conical flasks
- Pipettes
- Beakers
- Culture tubes
- Screw capped universal bottles
- Washing and cleaning – soap water, single distilled water, hot air oven dry, dust-proof cupboard/ drawer
- Cotton plug, aluminium foil
- Flame sterilization of instruments



2. Instruments

- Various sizes of scalpel
- Forceps
- Spatula
- Scissors, etc
- Should be of stainless steel

